

Better Detail For Better Quality
Better Detail For Better Quality



LINK CT-1000

Chair &

Stand Unit

New Technologies And Powerful Functions
Simple operation interface provides
with accurate operation experience.



◆ An automatic and full functional Chair & Stand Unit, automatic functions are realized in different parts: automatic table moving, automatic tableboard up and down, automatic column up and down, automatic phoropter arm.

◆ European smart design style in appearance, which looks simple and fashionable.

◆ Different instruments can be equipped on it: Auto Refractometer/Auto Phoropter/Slit Lamp/Auto Lensmeter/Trial Lens Set.

◆ Equipped with drawer for trial lens set, enough storage space is provided with.

◆ Wireless operation ways are provided with, such as Operation Panel/Mobile Control/Tablet Control.

◆ Leather chair with linear motor makes a quiet operation.

◆ An ideal Combi Unit for eye clinic, hospitals, optcial shops etc.

1. The table material is made of high-density film-coated board, which can be rotated by 90 degrees, which can place two instruments

2. The table can be moved left and right by electric and up and down by electric, and the operation is stable;

3. The column swing arm can move left and right electrically and up and down electrically, and runs smoothly;


4. Seat can be moved forward and backward, and electric reclining and electric lifting;

5. According to different customer groups, three different reset modes can be set;


6. Equipment color can be customized and change according to actual situation.



Tablet Control System: besides operation panel, it can be controlled by Tablet.



Automatic moving phoropter arm: left and right.



Aluminium Column moves up and down.



Automatic table moving: left and right, up and down.



Enough storage space for trial lens set etc.

Specification

| | |
|---|--|
| Maximum size: | 1200×1300×1800mm |
| Table size: | 950×450mm |
| Small table size: | 620×400mm |
| Drawer size: | 520×350mm |
| Control mode: | film button / phone remote control / flat panel remote control |
| Configuration equipment: | computer refractometer / electric cattle eye / slit lamp / focus meter |
| Table rotation angle: | 90° |
| Countertop left and right movement distance: | 400mm |
| Countertop up and down distance: | 120mm |
| Column (cattle eye arm) left and right moving distance: | 400mm |
| Column (cattle eye arm) up and down distance: | 120mm |
| Chair left and right movement distance: | 150mm |
| Chair up and down distance: | 150mm |
| Chair back angle: | 120-140 ° |
| Chair armrest load: | 100kg |
| Chair loading ability: | 300kg |
| Chair height: | 1400mm |
| Lamp: | LED (220 / 110V 9W) |
| Voltage: | AC220V / 110V 50 / 60HZ |
| Gross weight: | 220kg |
| Packing size: | 1600×1050×1020mm |
| Packaging volume: | 1.7 cubic |