



# LINK CT-500

## Chair and Stand Unit

New Technologies And Powerful Functions  
Simple operation interface provides  
with accurate operation experience.

# LINK CT-500

## Chair and Stand Unit

### LINK INSTRUMENTS

#### Features

- 1. Table material ABS hard plastic, can slide left and right, rotate 90 degrees outwards;
- 2. Place two instruments;
- 3. The seat has functions of rotation, aerodynamic reclining, movement, and electric lifting;
- 4. Mechanical swing arm feels good and accurate positioning;
- 5. The cabinet is small in size and stable in construction;
- 6. Rechargeable ophthalmoscope ophthalmoscope socket;
- 7. It can increase the matching use of the slit lamp, and charge according to the actual situation;
- 8. Seat color can be customized and charge according to actual situation;



Option1 The interface for slit lamp.



Option2 Don't need the interface.

#### Specification

Maximum size:	1100×800×2100mm
Table size:	890×405mm
Table moving angle:	90°
Table moving distance:	335mm
Swing arm angle:	+35°
Swing arm rotation range:	60°
Seat rotation angle:	45°
Seat height:	150mm (550–700mm)
Seat load:	250kg
Seat moving distance:	150mm
Seat armrest load:	100kg
Seat reclining angle:	120–140°
Light source:	LED (5V1A 3.6W)
Voltage:	AC220V / 110V 50 / 60HZ
Gross weight:	180kg
Packing size:	1200×1100×1050mm
Packing volume:	1.38 cubic

Better Detail For Better Quality